


If you prefer alginate impressions, may I suggest KROMOPAN 100. It is not time sensitive!

**KROMOPAN 100 – The original color-change alginate
Dust-free, Class Alginate Impression Material Gives you 100 hours of
dimensional stability Eliminates gagging, guess work and costly remakes A
perfect impression and a happier patient every time!**



Step one – Purple stage
Spatulate until pink color
appears.



Step two – Pink stage
Load tray and wait for
white color to appear.



Step three – White stage
Insert for impression (time in
mouth, 30 seconds).

Total working time: 1 minute 45 seconds from start to finish!